


# TITANIUM Total

CODE  
TAP 07.61


TI: Titanium ASTM F67 Gr. 2

Prosthesis with fixed length: from 3.00 mm to 8.00 mm

- WINDOWED (OPENED) FLANGE – it allows the clear view of the below surgery field.
- TITANIUM STRUCTURE – Titanium is a bioinert material, largely tested in surgery fields and with high mechanical resistance and lightness.



It would be better to put an autologous cartilage sheet between the prosthesis and neotympanum reducing possibility of encapsulation and rejection.



Wire  $\varnothing$ : 0.20 mm  
Columella  $\varnothing$ : 0.60 mm  
Weight: from 4,52 mg (TAP 07.61.300) to 5,22 mg (TAP 07.61.800)

Code	L (mm)
TAP 07.61.300	3.00
TAP 07.61.325	3.25
TAP 07.61.350	3.50
TAP 07.61.375	3.75
TAP 07.61.400	4.00
TAP 07.61.425	4.25
TAP 07.61.450	4.50
TAP 07.61.475	4.75
TAP 07.61.500	5.00
TAP 07.61.525	5.25
TAP 07.61.550	5.50
TAP 07.61.600	6.00
TAP 07.61.650	6.50
TAP 07.61.700	7.00
TAP 07.61.750	7.50
TAP 07.61.800	8.00

**Packaging**  
Dispenser box containing 6 sterile,  
individually packaged implants.