

Universal Adjustable Trimmable TITANIUM Prosthesis


CODE
TPL 07.57

Titanium GR. 2
PTFE: polytetrafluorethylene

This prosthesis consists of:

- ▶ titanium flange with 15 holes of 0.40 mm. diameter allowing penetration of connective tissue
- ▶ trimmable titanium stem
- ▶ universal PTFE columella for partial or total ossicular chain replacement.

In order to avoid extrusion, an autologous cartilage graft must be interposed between the flange of the prosthesis and the tympanic membrane.



Dimensions

Flange: mm 3.80 x mm 3.00

Total lengths

- of titanium stem: 8.00 mm
- of PTFE columella: 3.50 mm


Packaging

Dispenser box containing 6 sterile, individually packaged implants.

Code	Length	Columella Ø
TPL 07.57	Maximum length 8.50 mm (A)	0.70 mm
	Minimum length 3.50 mm (B)	1.80 mm (inner Ø 1.20 mm)

TPL 07.57 - Universal Adjustable Trimmable TITANIUM PROSTHESIS

Total version



TPL 07.57 - Universal Adjustable Trimmable TITANIUM PROSTHESIS

Partial version

